





MC-120A is designed ergonomically. It can be used with CO2 or N2O (the accessories are provided with the unit to use N2O). Cryo tips and the all probe are insulated so only the cryo tips frost. This unit works with high performance; the user obtains 3 mm depth in 30-45 sec. At the end of the operation; when the tumbler is left, thermal system works and defrosts in 4-5 sec. It is better to use it with the flexible hose. There is a button which is used to make the position of the tumbler stable on the cryo - probe. The user does not need to press the tumbler during the operation thanks to this system. There is an on-off valve on the cryo-probe so the cryo tips can be changed without turning off the CO2 / N2O cylinder. It has a carrier bag and a special cabinet which hides CO2 / N2O cylinder. It can be used for different purposes with a drawer, it also has trundles to be moved easily. There are a lot of special cryo tips (optional). MC-120A cryo unit has 2 years warranty for producing mistakes and 10 years warranty for the spares.





## **CRYOMED**

Tel:+90-312-278 19 46 Fax:+90-312-278 23 46 iletisim@bimedcryo.com www.bimedcryo.com





